



MITOCHONDRIAL ASSAYS

Mitochondrial Membrane Potential

TMRE/TMRM fluorescence responds to changes in ψ_m . TMRE fluorescence [excitation λ (Ex λ) of 568 nm and emission λ (Em λ) > 590 nm] in mitochondria will be visualized on a HCS (BD).

Inside a healthy, non-apoptotic cell, the lipophilic TMRE/TMRM dye, bearing a delocalized positive charge, enters the negatively charged mitochondria where it accumulates in an inner-membrane potential dependent manner. When the mitochondrial $\Delta\psi_m$ collapses in apoptotic cells, the TMRE /TMRM potentiometric dyes no longer accumulate inside the mitochondria and become more evenly distributed throughout the cytosol.

ROS

ROS generation will be detected with the fluorescence dye C-DCDHF-DA. The esterified fluorescence probe passes through cell membranes easily. Once inside the cell its lipophilic blocking groups are cleaved by non-specific esterases, resulting in a charged form that leaks out of the cells only slowly. The dye is useful for the detection of ROS. The plates were read in a FLIPRTETRA using excitation/emission (495 nm /520 nm)

Superoxides

Production of superoxide was evaluated intracellularly using the superoxide-sensitive dye DHE. DHE is oxidized by superoxide to a novel product which binds to DNA enhancing intracellular fluorescence. Cells were incubated with DHE (final concentration 10 μ M) in ensayo (NaCl 165 mM, KCl 4.5 mM, CaCl₂ 2 mM, MgCl₂ 1 mM, Hepes 10 mM y glucosa 10 mM, pH 7.4) for 30 min. The plates were read in a FLIPRTETRA using excitation/emission (515 nm /595 nm)

Mitochondrial Calcium

To determine mitochondrial calcium Rhod-2-AM will be used (Exc 549 / Em 578). Plates will be read on HCS (BD)

ATP/AMP/ADP levels

ATP, AMP and ADP concentrations will be determined by means of anion-exchange HPLC-UV using a binary gradient of 0.3 M ammonium carbonate and water on a ProPac PA1 column. ATP, AMP and ADP peaks will be detected compared their retention times.

NADH concentration

NADH oxidase activity will be measured as the aerobic conditions of 75 μ M NADH in the absence of external quinone substrates and other respiratory chain inhibitors.



Fundación MEDINA

Discovering the Future

Reactions rates will be calculated from the linear decrease of NADH concentration ($\lambda=340$ nm) in a spectrophotometer using submitochondrial particules.

Oxygen consumption

The assay is based on the ability of O₂ to quench the excited state of the MitoXpress probe. As the test material respire (e.g., isolated mitochondria, cell populations, small organisms, tissues and enzymes), O₂ is depleted in the surrounding solution/environment, which is seen as an increase in probe phosphorescence signal. Changes in oxygen consumption reflecting changes in mitochondrial activity are seen as changes in MitoXpress™ probe signal over time.

The assay is non-chemical and reversible, a decrease in oxygen consumption (an increase in O₂ levels) is seen as a decrease in probe signal. The assay will be read in a FLIPRTETRA

Fundación MEDINA

Parque Tecnológico de la Salud - Av. Del Conocimiento 3
18016 - Armilla - GRANADA
Tel. 958 993 965

www.medinaandalucia.es - info@medinaandalucia.es